

What it is

SampleVision™ v3.1 is a mobile friendly portal that allows any site, anywhere, to streamline the process of testing analytical samples. SampleVision connects remote sites and central labs seamlessly via Mirth Connect to any LIMS or LIS application on a smartphone, tablet, laptop, or PC, or will export requests to a spreadsheet. It's simple and intuitive to use, easy to set up, and fast to roll out.

Who it's for

LABS


- Analytical labs, that send out samples for analysis, staffed by non-laboratory personnel
- Central labs, that receive the sample test requests, and manage multiple remote labs

INDUSTRIES

- Chemical/Petrochemical
- Pharmaceutical Manufacturing
- Food and Beverage
- Clinical
- Environmental

How it works

SampleVision™ v3.1 securely submits sample information, prints barcodes, and documents chains of custody for samples to be sent to the testing laboratory. When testing is complete, lab results are sent via push notifications to access reports, with email sharing capabilities. Test requests can be sent through SampleVision via mobile app in remote sampling locations or through the web interface from **any device, anywhere, and at any time.**



Benefits

- ✓ Mobile friendly/cloud based
- ✓ Easily submit test orders from anywhere
- ✓ No LIMS license or knowledge needed
- ✓ Time savings
- ✓ Accurate and secure chain of custody
- ✓ Dashboard of test orders in the pipeline
- ✓ 24/7 access to results

Features

- Push notifications
- Part 11 compliance
- Customize/manage multiple accounts
- Automatic posting
- Multiple languages
- Geotagging capabilities
- Generates barcodes per sample
- Interfaces with any LIMS/LIS